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INSTRUCTIONS FOR USE - READ CAREFULLY!

Fresenius Kabi AG
61346 Bad Homburg
Germany

Purisole® SM diluted

Solution for Irrigation

Composition

1 liter contains:

Sorbitol	27.0 g
Mannitol	5.4 g
Water for injections in bulk, Ph. Eur.	
pH value	4.5 – 7.0
Titration acidity	< 1 mmol NaOH/l
Theoretical osmolarity	178 mosm/l

Indications

- Intra- and postoperative bladder irrigation in urological interventions such as transurethral resection of the prostate or bladder tumors.
- Distension and irrigation of the uterine cavity in operative hysteroscopy, e.g. endometrial ablation, septum dissection, resection of submucous myoma and synechiolysis.
- Distension and irrigation in operative arthroscopy.
- Distension in diagnostic procedures, e.g. diagnostic hysteroscopy, arthroscopy, and cystoscopy.

Contraindications

- Hypersensitivity to sorbitol and/or mannitol
 - Hereditary fructose intolerance: Purisole® SM diluted contains sorbitol. For this reason, the presence of hereditary fructose intolerance or a lack of fructose 1–6 diphosphatase must be
-
-

ruled out prior to irrigation with Purisole® SM diluted. Should anamnestic data be unavailable or unclear, a fructose intolerance test must be performed before irrigation. Non-observance of this contraindication may lead to potentially fatal complications.

- Hyperhydration
- Methanol intoxication
- Newborns and infants

Precautions

Caution is required in patients with pre-existing hyponatremia. Serum sodium levels must be corrected before use in this case.

Purisole® SM diluted should be used with caution in patients with cardiac or renal insufficiency.

Purisole® SM diluted should be used with caution in patients with diabetes mellitus (perioperative blood glucose monitoring recommended) or metabolic disorders (measuring of lactate levels recommended in the perioperative period).

Adequate precaution measures should be taken to avoid massive fluid absorption, e.g. balancing of the fluid in- and outflow (stopping the procedure if the difference exceeds 1 - 2 liters), limiting the operation time to 1 hour, and not exceeding the recommended maximum pressure for the intervention.

Side-effects

Fluid absorption may occur in operative trans-urethral, hysteroscopic, and arthroscopic procedures. In case of intravascular absorption of large amounts of electrolyte-free irrigation solutions, hyperhydration and changes in the electrolyte balance (hyponatremia) may occur and consequently lead to acid base imbalance (metabolic acidosis). This condition is known as „TUR syndrome“. The associated symptoms include nausea, vomiting, increase or decrease in blood pressure, bradycardia, headache, restlessness, confusion, lethargy, chest pain, muscular twitches, dyspnea,

and poor urine output. In severe cases seizures, shock, pulmonary edema, cerebral edema, cardiovascular collapse, and coma may occur. Fatal cases have been reported.

If the development of TUR syndrome is suspected, the procedure should be terminated, and adequate therapeutic measures initiated according to the individual condition of the patient.

Potential side effects associated with the metabolism of sorbitol are hyperglycemia in patients with diabetes mellitus and hyperlactatemia and/or lactic acidosis in patients with metabolic disorders or in case of massive absorption.

Interaction with other medicines

None known.

Dosage

Dosage and quantity to be applied will result from the type of intervention.

Description and connectivity

The Careflex® Duo bag is a polyolefin bag with a spiking port and the Care-Lock® port.

The Care-Lock® port is only compatible with the Care-Lock® transfer set.

The spiking port is compatible with all Fresenius Kabi transfer sets. Other spiking systems may be used and have to be tested for suitability by the user. Based on tests conducted with other transfer systems on the market the following sets are known to be not compatible with the spiking port of the Careflex® Duo bag: Pump set for ARTHREX-System by Meise Medizintechnik (Ref: AR 6410), Pump Tube for fms solo-duo by Meise Medizintechnik (Ref: AR 7410), Tubing for Cystoscope or Resectoscope by Polfa-Lublin, Normoflo Irrigation Warming Set by Smiths Medical (Ref: IRI-600), AHTO Disposable Suction Irrigator Tube Set with Tip by Stryker (Ref: 250-070-620), Flocontrol Tube Set Pump SIL, Reusable by Stryker (Ref: Z1428-39), TUR

Irrigation by Techno Med (Ref: TY30002SR), Fluido Irrigation URO Set by The Surgical Company International B.V. (Ref: 670800-B), Fluido Irrigation URO Spike Set by The Surgical Company International B.V. (Ref: 680803-B), Single Use Tubing Set for Pump PV-5201 by Vimex Sp.Z.o.o. (Ref: PV-5201SUTS).

Use of non-compatible transfer systems may lead to leakages or loss of spike during use and impossibility of irrigation.

An instruction how to spike the bag can be found at the end of this instruction for use.

Presentation

Polyolefin bags (Careflex® Duo) 6 x 1500 ml
4 x 3000 ml
2 x 5000 ml

Additional Information

Not for intravenous use or injection!

The sterility of Purisole® SM diluted is guaranteed only until the container is opened. The solution must be used immediately after opening.

Use only if the solution is clear and free from particles and if the package is undamaged.

For single use only.

Any unused portion of the solution must be discarded.

Re-use of medical single use equipment involves a risk of contamination of the equipment. This may cause infection in the patient or the user, which may lead to illness or death.


Keep out of the sight and reach of children!

Do not use Purisole® SM diluted after expiry date.

Date of revision

August 2021

	CE marking of conformity with Medical Devices Directive 93/42/EEC (as amended) with the intervention of Notified Body 0123
	Indicates the presence of a sterile fluid path with-in the medical device in cases when other parts of the medical device, including the exterior, might not be supplied sterile
	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Indicates the medical device manufacturer.
	Indicates the need for the user to consult the instructions for use.
	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Indicates the date after which the medical device is not to be used.
	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	Indicates a medical device that should not be used if the package has been damaged or opened.
	Indicates the manufacturer's catalogue number so that the medical device can be identified.

	This product is not for intravenous use or injection.
	Indicates a medical device that needs protection from light sources.
	Indicates a medical device that needs to be protected from moisture.
	Indicates the temperature limits to which the medical device can be safely exposed. Store at room temperature.


Spiking the bag


For spiking, the bag can be hung on an IV pole or alternatively placed on a flat surface. Administration of the content in hanging position..




Check product name, bag integrity and the expiry date. Inspect the container for damage or leakage and for particles. **Do not use if the container is damaged or if the solution contains particles.**

Identify the blue spiking port.

Break off the blue tamper-evident cover from the port.
Do not touch the opened surface!

Close clamp of transfer set if applicable. Insert the spike **centered** and **straight** until the blue plastic collar of the port meets the shoulder of the spike.

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